



D-TAIL LAB – TECHNICAL DATA SHEET (TDS)

MAVERICK — Premium Interior Multi-Purpose Cleaner — HI-Con Series

Document code: TDS-DMAV • Date: 2025-11-13

Quick Facts	
Product Name	MAVERICK Premium Interior Multi-Purpose Cleaner — HI-Con Series
Series	HI-Con
Color	Amber
Scent	Perfumed
pH (25°C)	7.0-8.0
Mobile Wash Friendly	Yes

Description & Key Benefits

Active-foaming, pro-strength interior cleaner for plastics, vinyl, coated leather, rubber, fabrics and carpets. Lifts grime and body oils for fast wipe-off; rinses clean and leaves a fresh, ready-to-dress finish. Balanced pH (7.0-8.0) for routine cabin maintenance.

Suitable / Do Not Use / Limitations

Suitable: dashboards, door cards, consoles, vinyl, coated leather, plastics, rubber, light fabrics/carpets.
Do Not Use: on hot surfaces or in direct sun; avoid saturating delicate leathers/alcantara; do not let dry.
Limitations: heavy stains may require a dedicated fabric cleaner (e.g., GAME OVER) or multiple passes.

Dilution Chart	
Light dirt (standard)	1:50
Average dirt	1:25
Strong cleaning	1:10

Directions for Use

Work on a cool interior, out of direct sun.
Spray onto the area or a brush/applicator; agitate gently.
Wipe residue with a clean microfiber; repeat as needed.
For delicate finishes, use a soft brush and avoid over-wetting.

Safety (Summary) — See SDS for full details	
GHS Word	Warning
Hazards	May cause skin and eye irritation. Harmful if swallowed.
Precautions	Wear gloves/eye protection. Avoid release to the environment. Keep out of reach of children.
First Aid	Eyes/Skin: Rinse with water. If irritation persists, get medical advice. If swallowed: rinse mouth and call a Poison Center.

Transport (Summary)	
UN / Proper Shipping Name	Not regulated as dangerous goods
Class / Packing Group	—
TDG / DOT / IMDG / IATA	Not regulated
Marine Pollutant	No

Environmental / Disposal

Avoid release to the environment.

Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage & Shelf Life

Store tightly closed in original container; protect from heat/freezing.

Shelf life: 24 months (unopened) when stored as directed.

Packaging & Item Numbers		
Size	SKU	Case Pack
500 mL	DMAV500	6
4 L	DMAV4LT	4
20 L	DMAV20L	

Manufacturer: D-TAIL LAB | Address: 200 Saint-Francois-Xavier, Unit 114D, Delson, QC | Phone: +1 514-931-4422 | Email: info@dtailab.com | Site: dtailab.com • See SDS for full details